

Aluminium powder		Al (26.98)
 DANGER	Flammable solid [H228]. In contact with water releases flammable gases [H261]. Aluminium powder is dangerous in contact with: • WATER. Under very hot conditions, hydrogen is formed. ACIDS and ALKALIS. Hydrogen is formed. • ALCOHOLS. Violent reactions occur after a lengthy induction period. • HALOGENS, ORGANOHALOGENS. SULFUR. Violent or explosive reactions occur. • METAL OXIDES, NITRATES(V), CHLORATES and other SOLID OXIDISING AGENTS. Explosive mixtures form. • MERCURY and its COMPOUNDS (includes solutions). Rapid oxidation of aluminium occurs, may cause ignition. WEL (mg m ⁻³): 10 (LTEL), 30 (STEL), as inhalable dust; 4 (LTEL), 12 (STEL), as respirable dust.	
	Storage Storage code: FW [Dull-grey powder] • An oxide layer forms on turnings and powder (which is particularly prone to oxidation). This surface oxide reduces the activity of the metal. Keep containers tightly closed and dry.	

Aluminium foil		Al (26.98)
Currently not classified as hazardous. Aluminium foil is dangerous in contact with: • ACIDS and ALKALIS. Hydrogen is formed. • HALOGENS. Violent reactions occur. • MERCURY and its COMPOUNDS (includes solutions). Rapid oxidation of aluminium occurs which may cause ignition		
Storage	Storage code: Gln	[Shiny-silver metal]

Emergencies	Follow standard procedures (see <i>Emergency Hazcards</i> and GL 120). NOTE: • If aluminium powder catches fire: Although difficult to ignite, the powder burns vigorously and is difficult to extinguish. If safe to do so, extinguish small fires by smothering with clean, dry sand (keep a container of sand specifically for this purpose).
--------------------	---

Detailed guidance on specific activities and techniques involving these substances can be found at: www.cleapss.org.uk

General use of:	Hazard information	User	Suggested general control measures and guidance
Aluminium powder	 DANGER (See reverse)	TT (Y7)	- Wear eye protection. - Avoid raising dust. Keep powder away from naked flames. Other notes; see also detailed activity-specific guidance: - Thermite reaction (TT): See <i>Thermite reaction</i>. - Powdered aluminium is particularly dangerous if blown into a Bunsen burner flame. Do not use aluminium powder in reactions with acids, alkalis, chlorine or bromine; use small pieces of foil. Disposal: W2, WSpec; see below
Aluminium foil	Currently not classified as hazardous but see notes on reverse.	Y7	Eye protection and/or other control measures may be advised for some procedures. Check activity-specific guidance and/or GL 120. Disposal: W3, W8
Disposal	Follow general guidance in <i>About Hazcards</i> (GL 120). NOTE for <u>small amounts</u> of aluminium powder: WSpec: Work in a fume cupboard. Ensure no naked flames. Add 10 g of the powder to 1 litre of 2 M hydrochloric acid (in a bucket or large beaker). Ensure solution remains acidic (use an indicator); add more acid if necessary. Leave mixture in the fume cupboard for an hour or until reaction appears to be complete. Pour resulting solution down a foul-water drain, diluting it further with a bucket of water.		